

WORKSTATION CAMBRIDGE


Timeless, adaptable & *versatile*.

.....





Scan the QR code to see this product on our website or [click here](#)



W (width, mm)	D (depth, mm)	h (worktop height, mm)	H (column height, mm)
1.500	750	900	1.570
1.200	900	750	1.800
900			

MAX Load 350 Kg

Power supply side can be switched left to right and viceversa. Check [accessories pages](#) to see how shelves can be added, removed and changed position.

- Underbench shelf available.
- Gas taps can be attached to column.
- Some other options available.
- In addition, special accessory attachments can be design and manufactured on demand.

Main frame fully welded steel, epoxy coated 9010 (check other options) with M5 nuts below worktop for accesories and attachments.

Columns manufactured from metal sheet, magnetic covers, all joints rubber-sealed. Top covers removable allows for different combinations. Epoxy coated 9010 (check other options) with M8 nuts at internal faces for accesories. Nuts in front faces as an option.

Worktop material on demand. Advised HPL 16 mm high chemical resistance.

Feet, articulated with rubber base, interchangeable with combi castor-foot.

Data setup on demand, advised to be magnet clip-routed.

Electrical setup 8 power sockets installed (3 each side above worktop). **Circuit breaker** setup on demand.

