

EQUIPMENT BENCH

Engineered for
strength, stability, and confidence.

.....




Equipment Bench

Engineered for *strength, stability, and confidence.*



Scan the QR code to see this product on our website or [click here](#)



W (width, mm)	D (depth,mm)	h (worktop height, mm)
1.500	600 / 750 / 900	900 / 750
1.200	600 / 750 / 900	900 / 750
900	600 / 750 / 900	900 / 750
750	600 / 750 / 900	900 / 750

Max recommended load - 400 kg


Check **accessories** pages to see how:

- **Underbench shelf** available.
- Some **other options** available.
- In addition, **special accessory attachments** can be design and manufactured on demand.

Main frame fully welded steel, epoxy coated 9010 (check other options) **with M5 nuts below worktop** for accesories and attachments.

Worktop material on demand. Advised HPL 16 mm high chemical resistance.

Feet, articulated with rubber base, interchangeable with combi castor-foot.



Equipment Bench is ideal as the head of an island setup. *(Highlighted in green)*



As a **back to back position with one of our WorkStation**, you can get a simpler yet fully working setup, with a wider working surface (perfect for oversize equipment). *(Highlighted in green)*



750x750mm Equipment Bench is specially designed to solve 90° WorkStation to WorkStation setup, without cornering these, keeping mobility agile. *(Highlighted in green)*



Technical Data Sheets

Visit our
website for
more info



www.ibr-innovasystem.com

Contact us at: hello@ibr-innovasystem.com

