

PRIMARY BENCH 4P

*A clean, adaptable, and efficient
solution for everyday lab work.*

.....




Primary Bench 4P

A clean, adaptable, and efficient solution for everyday lab work.



Scan the QR code to see this product on our website or [click here](#)



W (width, mm)	D (depth, mm)	h (worktop height, mm)
1.500	600 / 750 / 900	900 / 750
1.200	600 / 750 / 900	900 / 750
900	600 / 750 / 900	900 / 750
750	600 / 750 / 900	900 / 750

Max recommended load - 325 kg

Accessories and attachments such as electrical trunking brackets, side and back covers, **or on demand new designs, can be attached** at back beam and back legs nuts.


Primary Bench can be configured with appliances, check Primary+ Bench sheet for more info.

Check Primary Bench 5P sheet for info on **extralarge sizes (up to 2.400 mm width)**.

Main frame fully welded steel, epoxy coated 9010 (check other options) **with M5 nuts below worktop (back beam) and M6 at back legs** for accesories and attachments.

Worktop material on demand. Advised HPL 16 mm high chemical resistance.

Feet, articulated with rubber base.




Primary Benches can also **form islands and peninsulas on their own**. (Highlighted in green).



Primary Bench is an ideal solution for a traditional approach in **fixed positions, small rooms, corners, etc.**



Makes a **great complement for a comb layout** setting the wall fixed benching as a starting point **for comb's peninsulas, marking widths and spacing perfectly**, by combining 1.500 mm width Primary Benching with 750 mm deep WorkStations. (Highlighted in green).



Technical Data Sheets

Visit our
website for
more info



www.ibr-innovasystem.com

Contact us at: hello@ibr-innovasystem.com

