

PRIMARY BENCH 5P

All welded, lightweight
frame, *extra leg*.

.....




Primary Bench 5P

Wide surface to work,
power to perform.



Scan the QR code
to see this product
on our website or
[click here](#) 



W (width, mm)	D (depth, mm)	h (worktop height, mm)
1.800	600 / 750 / 900	900 / 750
2.100	600 / 750 / 900	900 / 750
2.400	600 / 750 / 900	900 / 750

Max recommended load - 325 kg


Accessories and attachments such as electrical trunking brackets, side and back covers, **or on demand new designs, can be attached** at back beam and back legs nuts.

Primary Bench can be configured with appliances, check Primary+ Bench sheet for more info.

Main frame fully welded steel, epoxy coated 9010 (check other options) **with M5 nuts below worktop** (back beam) **and M6 at back legs** for accesories and attachments.

Worktop material on demand. Advised HPL 16 mm high chemical resistance.

Feet, articulated with rubber base.



Primary Benches can also conform islands and peninsulas on their own. Wider islands can be arranged with Primary Bench 5P, having a bigger worktop area for lighter tasks that requires more space.




Primary Bench goes perfect for a **traditional approach in fixed positions, small rooms, corners, etc.** Version with five posts adds extra leg room.



Technical Data Sheets

Visit our
website for
more info



www.ibr-innovasystem.com

Contact us at: hello@ibr-innovasystem.com

