


CABINET ON CASTORS

A reliable storage essential,
designed to complement your workflow.

.....




Cabinet On Castors

A reliable storage essential, *designed to complement your workflow.*



Scan the QR code to see this product on our website or [click here](#) 



W (width, mm)	D (depth, mm)	h (worktop height, mm)
450	550	830 / 660
600	550	830 / 660
850	550	830 / 660
1.200	550	830 / 660

MAX Load 25 kg - per drawer
MAX Load 100 kg - total

.....

Can be delivered in **fully HPL or Melamine.**

Melamine version **also has HPL worktop.**

Optional **270° opening hinges.**

Check availability for **other colours or door-drawer combinations** out of this list.

Body manufactured fully with **16 mm melamine with black PVC edging (advised)** or 16 mm HPL.

Usable worktop built on 16mm HPL usable top.


White finishing as standard (check options).

Rubber wrapped plastic core castors, two without brakes at back, two with brakes at front.


Black handles (200 mm).

Soft closing hinges opening 165°. Optional 270° opening hinge.

Models with door(s) includes **height adjustable internal shelf.**



Underbench configuration examples



Possible combinations



Available **door and drawers combinations**
for **single door** cabinets (widths 450 and 600 mm)



Available **door and drawers combinations**
for **double door** cabinets (widths 850 and 1.200 mm)

Technical Data Sheets

Visit our
website for
more info



www.ibr-innovasystem.com

Contact us at: hello@ibr-innovasystem.com

