

WORKSTATION BARCELONA

Our most iconic Workstation;
Versatile. Distinctive. Refined.


.....



Our most iconic WorkStation;
Versatile. Distinctive. Refined.



Scan the QR code
to see this product
on our website or
[click here](#)



W (width, mm)	D (depth, mm)	h (worktop height, mm)	H (column height, mm)
1.500	750 / 900	750 / 650	1.570 / 1.800
1.200	750 / 900	750 / 650	1.570 / 1.800
900	750 / 900	750 / 650	1.570 / 1.800
750	750 / 900	750 / 650	1.570 / 1.800

Max recommended load - 400 kg

Power supply side can be switched left to right and viceversa. Check **accessories pages to see how** shelves can be added, removed and changed position.

- **Underbench shelf** available.
- **Gas taps** can be attached to column.
- Some **other options** available.
- In addition, **special accessory attachments** can be design and manufactured on demand.

Main frame fully welded steel, epoxy coated 9010 (check other options) **with M5 nuts below worktop** for accesories and attachments.


Columns manufactured in **steel tube bolted to main frame, epoxy coated 9010** (check other options), with **M8 nuts at internal faces** for accesories. **Nuts in front faces** as an option.

Worktop material on demand. Advised HPL 16 mm high chemical resistance.

Feet, articulated with rubber base, interchangeable with combi castor-foot.


Data setup on demand, advised to be magnet clip-routed.

Electrical setup: 8 power sockets installed (4 each side above worktop). **Circuit breaker setup** on demand.



Technical Data Sheets

Visit our
website for
more info



www.ibr-innovasystem.com

Contact us at: hello@ibr-innovasystem.com

