

ATB BENCH


All Time Bench.

.....





Scan the QR code
to see this product
on our website or
[click here](#)



All dimensions available on request



Examples:

W (width, mm)	D (depth, mm)	h (worktop height, mm)
1.500	750	900
1.000	750	450
3.000	750	900

Max recommended load - 325 kg

Ask for special setups such as:

- Underbench shelf.
- Undermounted full cabinet.
- Extra legs / Reinforcements.

Main frame fully welded steel, epoxy coated 9010 or 9016 (check other options).

Worktop material available on request.
Advised HPL 16 mm high chemical resistance.

Feet:

- Articulated with rubber base /
- Industrial-grade castors

On demand. Not interchangeable.



Main frame manufactured in **square section steel tube**.



Some examples of **worktop materials**.

CONFIGURABLE LOWER BEAM

Lower beam position on demand:

- **Back**, for wall position.
- **Middle**, for island or mobile position.



Leveling Feet or Industrial Grade Castor
On demand. Not interchangeable.



Leveling feet



Industrial-grade castor

Made-to-Measure

Our **highly scalable architecture** accommodates any requirement, bridging the gap between a **450mm micro-station** and a **heavy-duty 3-meter configuration**.



Extra large work surface, with middle lower beam and Industrial-grade castors.



Standard Mobile Unit with equipment storage.



Customised Bench:

- Underbench reinforced shelf
- 1 door cabinet

Micro-Station with middle lower beam.



Technical Data Sheets

Visit our
website for
more info



www.ibr-innovasystem.com

Contact us at: hello@ibr-innovasystem.com

