

NOVA AIR

Ductless Fume Cupboard

*A breath of fresh air
in laboratory safety.*


.....



Ductless Fume Cupboard
+ Equipment Bench



Scan the QR code
to see this product
on our website or
[click here](#)




Mechanical Data Airone R3 Ducted Fume Cabinet	
Min. Working Height	200 mm
Max. Working Height	365 mm
Usable Internal Height	780 mm
Power Supply	230V, 50 Hz
Main Construction	Mild Steel
Stand	Mild Steel

W (width, mm)	D (depth, mm)	H (column height, mm)	h (worktop height, mm)	INTERNAL (WxDxH, mm)
1.000	760	2.200	900	950x550x780
1.200	760	2.200	900	1.150x550x780
1.500	760	2.200	900	1.450x550x780
1.800	760	2.200	900	1.750x550x780

- Recirculating Technology** eliminates the need for external ductwork.
- Digital Control System** provides real-time monitoring of airflow, filter status, and system performance.
- Carbon / HEPA Filter configurations**, customisable to suit specific chemical applications.
- Low Energy Digital Fans** that deliver consistent airflow while reducing power consumption.
- Low Fan Noise** creates a **quieter working environment** that supports concentration and comfort.

- Large Front Aperture** allows for **easy access to the working area**, facilitating smooth handling of equipment and materials.
- Low Airflow Alarm** provides an **immediate alert** if airflow drops below safe levels, helping **prevent accidental exposure** to hazardous fumes or chemicals.
- Hinged Front Panel** aids routine maintenance, allowing **quick access to filters and internal components**.
- Integral LED Lighting** is **integrated into the unit**, providing clear visibility across the working area without adding heat or glare.
- Filter saturation Alarm**, alert users when filters are approaching capacity.



Equipment Bench

Integrated bench solution, which offers flexibility in lab layout and mobility.



Vented Storage Cabinets provide secure, ventilated storage for chemicals and reagents, helping to **reduce the risk of vapour build-up and improve overall lab safety.**

Special Bespoke Units designed to your exact specifications: dimensions, specialist filtration, or integrated utilities.

Technical Data Sheets

Visit our
website for
more info



www.ibr-innovasystem.com

Contact us at: hello@ibr-innovasystem.com

